

PDD-RELATED MEDICAL DEVICES

PMA-APPROVED REPROCESSING METHODS

Reprocessing Methods			Rigid Telescope System			Flexible Endoscope System
			Tricam® PDD Camera Heads	Rigid HOPKINS® Telescopes	Fluid Light Cable	Flexible PDD Video Cystoscopes
			20221139/20221137	27005 AIA/BIA/GIA/FIA	495 FS	11272 VNIUA/11272 VNIA
Sterilization	Steam	Pre-Vacuum, 132°C/4 min.	N/A	✓ Approved	N/A	✓ Approved - PDD Video Cystoscope Suction Unit Assembly 11272 AE
	Ethylene Oxide	100% EO	✓ Approved	✓ Approved	N/A	Not needed
	STERRAD® 100S		✓ Approved	✓ Approved	N/A	N/A
	STERRAD® NX®	Standard Cycle	✓ Approved*	✓ Approved		N/A
		Advanced Cycle	N/A	N/A	✓ Approved	✓ Approved**
	STERRAD® 100NX®	Standard Cycle	✓ Approved*		N/A	N/A
		Flex Cycle	N/A	N/A	N/A	✓ Approved**
		DUO Cycle	✓ Approved*	N/A		✓ Approved**
	V-PRO® maX	Non-Lumen Cycle				N/A
		Lumen Cycle				N/A
		Flexible Cycle	N/A	N/A	N/A	✓ Approved**
	V-PRO® 60	Non-Lumen Cycle	✓ Approved*			N/A
		Lumen Cycle	✓ Approved*			N/A
		Flexible Cycle	N/A	N/A	N/A	✓ Approved**
STERIS® SYSTEM 1E®			N/A		✓ Approved*	
High-Level Disinfection	CIDEX®	2.4% Glutaraldehyde Soln • 45 minutes	✓ Approved	✓ Approved	✓ Approved	
	Revital-Ox™ RESERT® HLD	2% Hydrogen Peroxide Soln • 8 minutes				✓ Approved*

* On 10/23/2018

** On 2/21/2019

Revital-Ox is a trademark and RESERT, STERIS, V-PRO, SYSTEM 1E, and SYSTEM 1 are registered trademarks of STERIS Corporation.

CIDEX, STERRAD, NX and 100NX are registered trademarks of Cilag GmbH International.



800-421-0837

www.karlstorz.com